

Table 6
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Method	Laboratory Reporting Limits ²	Holding Time
Non-Preserved Inorganics - 2-1 Liter Plastic				
Alkalinity	Cool < 4°C	SM 2320B	5.0 mg/L	14 Days
Chloride	Cool < 4°C	SW 9056	1.0 mg/L	28 Days
Sulfate	Cool < 4°C	SW 9056	1.0 mg/L	28 Days
Biochemical Oxygen Demand (BOD)	Cool < 4°C	SM 5210B	3.0 mg/L	48 Hours
Preserved Inorganics - 1-500 mL Plastic				
Nitrate+Nitrite	H ₂ SO ₄ , Cool < 4°C	SM 4500 NO3 E	0.1 mg/L	28 Days
Ammonia	H ₂ SO ₄ , Cool < 4°C	EPA 350.1	0.2 mg/L	28 Days
Total Organic Carbon (TOC)	H ₂ SO ₄ , Cool < 4°C	SM 5310C	1.0 mg/L	28 Days
Chemical Oxygen Demand (COD)	H ₂ SO ₄ , Cool < 4°C	SM 5220D	5.0 mg/L	28 Days
Preserved Inorganics - 1-500 mL Plastic				
Sulfide	NaOH & Zn Acetate, Cool < 4°C	SM 4500 S2 C,F	2.0 mg/L	7 Days
Total Metals - 1-500 mL Plastic				
Iron	HNO ₃ , Cool < 4°C	SW 6010C	0.1 mg/L	6 Months
Manganese	HNO ₃ , Cool < 4°C	SW 6010C	0.01 mg/L	6 Months
Calcium	HNO ₃ , Cool < 4°C	SW 6010C	1.0 mg/L	6 Months
Potassium	HNO ₃ , Cool < 4°C	SW 6010C	1.0 mg/L	6 Months
Sodium	HNO ₃ , Cool < 4°C	SW 6010C	1.0 mg/L	6 Months
Dissolved Metals - 1-250 mL Plastic				
Iron	HNO ₃ , Cool < 4°C	SW 6010C	0.1 mg/L	6 Months
Manganese	HNO ₃ , Cool < 4°C	SW 6010C	0.01 mg/L	6 Months
RCRA App I VOCs - 3-40 mL VOA Vials				
Acetone	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Acrylonitrile	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Benzene	HCl, Cool < 4°C	SW 8260B	0.5 ug/L	14 Days
Bromochloromethane	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Bromodichloromethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Bromoform	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Bromomethane	HCl, Cool < 4°C	SW 8260B	4.0 ug/L	14 Days
2-Butanone (MEK)	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Carbon disulfide	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Carbon tetrachloride	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
Chlorobenzene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Chlorodibromomethane	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Chloroethane	HCl, Cool < 4°C	SW 8260B	4.0 ug/L	14 Days
Chloroform	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Chloromethane	HCl, Cool < 4°C	SW 8260B	3.0 ug/L	14 Days
1,2-Dibromo-3-chloropropane	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
1,2-Dibromoethane	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Dibromomethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,2-Dichlorobenzene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,4-Dichlorobenzene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1-Dichloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,2-Dichloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1-Dichloroethene	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
cis-1,2-Dichloroethene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
trans-1,2-Dichloroethene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
trans-1,4-Dichloro-2-butene	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days

Table 6 (cont)
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Methods	Laboratory Reporting Limits ²	Holding Time
RCRA App I VOCs - 3-40 mL VOA Vials				
1,2-Dichloropropane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
cis-1,3-Dichloropropene	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
trans-1,3-Dichloropropene	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Ethylbenzene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
2-Hexanone	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Iodomethane	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
4-Methyl-2-pentanone (MIBK)	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Methylene chloride	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Styrene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1,1,2-Tetrachloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1,1,2,2-Tetrachloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Tetrachloroethene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Toluene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1,1-Trichloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,1,2-Trichloroethane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Trichloroethene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Trichlorofluoromethane	HCl, Cool < 4°C	SW 8260B	4.0 ug/L	14 Days
1,2,3-Trichloropropane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Vinyl Acetate	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
Vinyl Chloride	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
m+p-Xylenes	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
o-Xylene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
RCRA App II Inorganics - 1-500 mL Plastic				
Cyanide	NaOH, Cool < 4°C	SW 9012A	0.01 mg/L	14 Days
RCRA App II VOCs - 3-40mL VOA Vials				
Acetonitrile	Cool < 4°C	SW 8015m	10.0 mg/L	14 Days
Acrolein	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days
Allyl chloride	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
Chloroprene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
1,3-Dichlorobenzene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Dichlorodifluoromethane	HCl, Cool < 4°C	SW 8260B	3.0 ug/L	14 Days
1,3-Dichloropropane	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
2,2-Dichloropropane	HCl, Cool < 4°C	SW 8260B	4.0 ug/L	14 Days
1,1-Dichloropropene	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Ethyl methacrylate	HCl, Cool < 4°C	SW 8260B	2.0 ug/L	14 Days
Hexachlorobutadiene	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Isobutanol	Cool < 4°C	SW 8015	10.0 ug/L	14 Days
Methacrylonitrile	HCl, Cool < 4°C	SW 8260B	1.0 ug/L	14 Days
Methyl methacrylate	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Naphthalene	HCl, Cool < 4°C	SW 8260B	5.0 ug/L	14 Days
Propionitrile	HCl, Cool < 4°C	SW 8260B	10.0 ug/L	14 Days

Table 6 (cont)
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Methods	Laboratory Reporting Limits ²	Holding Time
RCRA App II SVOCs - 2-1 Liter Glass				
Acenaphthene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Acenaphthylene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Acetophenone	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Acetylaminofluorene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Aminobiphenyl	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Anthracene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzo(a)anthracene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzo(b)fluoranthene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzo(k)fluoranthene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzo(g,h,i)perylene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzo(a)pyrene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Benzyl alcohol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
bis(-2-chloroethoxy)Methane	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
bis(-2-chloroethyl)Ether	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
bis(2-chloroisopropyl)Ether	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
bis(2-ethylhexyl)Phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Bromophenyl phenyl ether	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Butyl benzyl phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Chloroaniline (p-Chloroaniline)	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Chlorobenzilate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Chloro-3-methylphenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Chloronaphthalene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Chlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Chlorophenyl phenyl ether	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Chrysene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
3/4-Methylphenol (m+p-Cresols)	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Methylphenol (o-Cresol)	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Diallate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Dibenzo(a,h)anthracene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Dibenzofuran	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Di-n-butyl phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
3,3'-Dichlorobenzidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,4-Dichlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,6-Dichlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Diethyl phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Thionazin	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
7,12-Dimethylbenz(a)anthracene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Dimethoate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
p-(dimethylamino)Azobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
3,3'-Dimethylbenzidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,4-Dimethylphenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Dimethyl phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,3-Dinitrobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4,6-Dinitro-2-methylphenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,4-Dinitrophenol	Cool < 4°C	SW 8270C	20.0 ug/L	7 Days
2,4-Dinitrotoluene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,6-Dinitrotoluene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days

Table 6 (cont)
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Methods	Laboratory Reporting Limits ²	Holding Time
RCRA App II SVOCs - 2-1 Liter Glass				
Di-n-octyl phthalate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Diphenylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Disulfoton	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Ethyl methanesulfonate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Famphur	Cool < 4°C	SW 8270C	20.0 ug/L	7 Days
Fluoranthene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Fluorene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Hexachlorobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Hexachlorocyclopentadiene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Hexachloroethane	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Hexachloropropene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Indeno(1,2,3-cd)pyrene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Isodrin	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Isophorone	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Isosafrole	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Kepone	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Methapyrilene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
3-Methylcholanthrene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Methyl methanesulfonate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Methylnaphthalene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Methyl parathion	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,4-Naphthoquinone	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1-Naphthylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Naphthylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Nitroaniline	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
3-Nitroaniline	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Nitroaniline	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Nitrobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2-Nitrophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
4-Nitrophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitroso-di-n-butylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosodiethylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosodimethylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosodiphenylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitroso-di-n-propylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosomethylethylamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosopiperidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
n-Nitrosopyrrolidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
5-Nitro-o-toluidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Parathion	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Pentachlorobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Pentachloronitrobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Pentachlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Phenacetin	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Phenanthrene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Phenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,4-Phenylenediamine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days

Table 6 (cont)
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Methods	Laboratory Reporting Limits ²	Holding Time
RCRA App II SVOCs - 2-1 Liter Glass				
Phorate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Pronamide	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Pyrene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Safrole	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,2,4,5-Tetrachlorobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,3,4,6-Tetrachlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
o-Toluidine	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,2,4-Trichlorobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,4,5-Trichlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
2,4,6-Trichlorophenol	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
o,o,o-Triethyl phosphorothioate	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
1,3,5-Trinitrobenzene	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
Dinoseb	Cool < 4°C	SW 8270C	10.0 ug/L	7 Days
RCRA App II Pesticides - 2-1 Liter Glass				
Aldrin	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
alpha-BHC	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
beta-BHC	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
delta-BHC	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
gamma-BHC (Lindane)	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
alpha-Chlordane	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
gamma-Chlordane	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Chlordane	Cool < 4°C	SW 8081A	2.00 ug/L	7 Days
Dieldrin	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
4,4'-DDD	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
4,4'-DDE	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
4,4'-DDT	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Endosulfan I	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Endosulfan II	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Endosulfan sulfate	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Endrin	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Endrin aldehyde	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Heptachlor	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Heptachlor epoxide	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Methoxychlor	Cool < 4°C	SW 8081A	0.032 ug/L	7 Days
Toxaphene	Cool < 4°C	SW 8081A	2.0 ug/L	7 Days

Table 6 (cont)
Semiannual Groundwater Monitoring Parameters
Sioux Falls Regional Sanitary Landfill

Parameters & Sample Container ¹	Preservatives	Methods	Laboratory Reporting Limits ²	Holding Time
RCRA App II PCBs - 1-1 Liter Glass				
PCB 1016	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1221	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1232	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1242	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1248	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1254	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
PCB 1260	Cool < 4°C	SW 8082	0.8 ug/L	7 Days
RCRA App II Herbicides - 2-1 Liter Glass				
2,4-D	Cool < 4°C	SW 8151A	1.0 ug/L	7 Days
2,4,5-T	Cool < 4°C	SW 8151A	1.0 ug/L	7 Days
2,4,5-TP (Silvex)	Cool < 4°C	SW 8151A	1.0 ug/L	7 Days

Table adapted from Sioux Falls Regional Sanitary Landfill Comprehensive Monitoring Plan, September and November 2011

mg/L - milligrams per liter

ug/L - micrograms per liter

S.U. - standard units

¹Container listed is typical. Container selected for sample collection to be approved by laboratory performing analysis for the parameter listed.

²Laboratory reporting limits with no dilutions